ACRYLIC RESIN FOR PATTERNS

- (ES) Resina acrílica para modelos
- (FR) Résine pour modèles
- (DE) HS-Modellierkunststoff
- (IT) Resina acrilica calcinabile per modelli
- (NL) Modelleerkunststof
- (PT) Resina acrílica para modelos
- (CS) Akrylová modelovací pryskyřice

INSTRUCTIONS FOR USE

Areas of application

Manufacture and modelling of

- telescopic and conus crowns
- individual attachments
- fixations / temporary splints for soldering
- inlays / onlays / connectors and clasps
- Maryland bridges
- direct intraoral core build-ups
- implant superstructures

Contraindications:

In the event of proven allergy to one or more of the product ingredients, the product should not be used. If the product has come into direct contact with the oral mucosa, rinse the mouth thoroughly with water.

Instructions for use – Telescopic and Conus crowns

Preparation: Place the correct amount of water and liquid in separate beakers.

Processing: Dip the tip of the brush into the liquid and then pick up a small amount of powder. This can then be applied directly to the model and repeated, progressively building the secondary crown. Curing time: approximately 3 minutes depending on ambient temperature.

After completion of the secondary crown, carefully remove and check the margin and inside surface. Reposition this on the primary crown and carefully grind to a uniform thickness of approximately 0.3 - 0.4mm with an suitable bur. Embedding, burning out and casting are carried out in the same way as the waxing technique.

Advice:

- If a plaster model is being worked on, isolation should be used as separating agent
- Clean the brush with the liquid either during with work procedure or after use
- Do not return leftover liquid or powder to the bottle

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HENRY SCHEIN SERVICES GmbH